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
[https://www.nms.ac.jp/college/pickup\\_contents/pickup25\\_ishi.html](https://www.nms.ac.jp/college/pickup_contents/pickup25_ishi.html)

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<https://www.ppmx.com/business/products/>

## Supplementary Material

### Optimized protocols for immunohistochemical staining of rat and mouse sections using PPZ0506 antibody, with some modifications

#### Tissue preparation

1. Euthanize animals by intraperitoneal injection of mixed anesthesia (medetomidine, midazolam, and butorphanol).
2. Perfuse transcardially with 0.9% (w/v) saline and 4% (w/v) paraformaldehyde in 0.1 M phosphate buffer (PB; pH, 7.4).
3. Post-fix organs with the same fixative overnight at 4°C.

#### Preparation of paraffin-embedded sections of peripheral organs

1. Immerse fixed organs in PB and dehydrate in an ethanol series (70%, 80%, 90%, 95%, 99%, and 100%) and xylene.
2. Embed the fixed organs in paraffin.
3. Slice paraffin-embedded organs at 5- $\mu$ m thickness on a microtome (RM2235; Leica Biosystems, Nußloch, Germany).
4. Mount sections on Frontier-coated glass slides (Matsunami Glass, Osaka, Japan).
5. Deparaffinize with xylene and hydrate in an ethanol series (100%, 99%, 95%, 90%, 80%, and 70%).
6. Rinse with distilled water (DW) for 5 min.

#### Preparation of frozen sections of brains

1. Immerse fixed organs in 0.1 M PB containing 20% (w/v) sucrose for 2–3 days at 4°C.
2. Freeze the fixed organs in n-hexane at -80°C.
3. Slice frozen organs at a thickness of 30  $\mu$ m (for mouse brain) or 40  $\mu$ m (for rat brain) on a cryostat (CM3050S; Leica Biosystems).
4. Rinse three times with 0.1 M phosphate-buffered saline (PBS; pH, 7.4) for 5 min each time.
5. Mount sections on Frontier-coated glass slides.
6. Dry the slides on a slide warmer (PS-53; Sakura Finetek Japan, Tokyo, Japan) at 60°C for 1 h.
7. Rinse with DW for 5 min.

#### Immunohistochemical staining of rat and mouse sections with PPZ0506

1. Place sections mounted on glass slides in methanol containing 0.3% (v/v) H<sub>2</sub>O for 5 min.
2. Rinse with DW and PBS for 5 min each.

3. Autoclave at 121°C for 10 min in citrate-based antigen unmasking solution (H 3300; Vector Laboratories, Newark, CA, USA) for antigen retrieval.
4. Rinse three times with PBS (5 min for each rinse).
5. Blocking Rat: Block with 5% normal goat serum (NGS) in 0.1 M PBS containing 0.1% Triton X-100 (PBST; pH, 7.4) for 30 min. Mouse: Block with blocking reagent (component of Mouse-on-Mouse Polymer IHC Kit; ab269452; Abcam, Cambridge, UK) for 30 min.
6. Primary antibody reaction\* Rat: React with PPZ0506 antibody (1:2000 for brain tissues, 1:4000 for peripheral tissues; Perseus Proteomics, Tokyo, Japan) in PBST containing 5% NGS at 4°C overnight. Mouse: React with PPZ0506 antibody (1:2000 for brain tissues, 1:4000 for peripheral tissues) in PBST at 4°C overnight. \*To obtain negative controls, prepare sections subjected to solvent without primary antibody.
7. Rinse three times with PBST (5 min for each rinse).
8. Secondary antibody reaction Rat: React with 50% (v/v) goat anti-mouse immunoglobulin G conjugated to a horseradish peroxidase (HRP)-labeled polymer (ab214879; Abcam) in PBST for 2 h. Mouse: React with HRP polymer detector reagent (component of Mouse-on-Mouse Polymer IHC Kit; ab269452; Abcam) for 15 min.
9. Rinse three times with PBS (5 min for each rinse).
10. Rinse with 0.05 M Tris-HCl (pH, 7.5) for 5 min.
11. Visualize immunoreactive signals with 200 µg/mL of 3,3'-diaminobenzidine tetrahydrochloride (Merck KGaA, Darmstadt, Germany) and 0.01% (v/v) H<sub>2</sub>O<sub>2</sub> in 0.05 M Tris-HCl (pH, 7.5).
12. Counterstain with hematoxylin for peripheral tissues.
13. Rinse with DW, dehydrate in ethanol series (70%, 80%, 90%, 95%, 99%, and 100%), and clear in xylene.
14. Coverslip with Permount Mounting Medium (Thermo Fisher Scientific, Waltham, MA, USA)